

# PERALEX 9 Hecto +

PERACETIC ACID-BASED  
READY-TO-USE DISINFECTANT



## ■ APPLICATION

- Operating room, endoscopy, optical equipment

## ■ PRODUCT BENEFITS

- Ready-to-use (without activator)
- High-level disinfection in 10 minutes
- Compatibility with materials used in endoscope disinfection benches, in particular corian
- Formulated without aldehydes: no fixation of proteins



**Fb Franklab®**  
notre expertise l'Ultra-Propreté



READY  
TO USE

CE  
0459

# PERALEX 9 Hecto +

## PERACETIC ACID-BASED READY-TO-USE DISINFECTANT

### INSTRUCTION FOR USE

#### Preparation of the disinfection bath :

- Wear reusable long-sleeved gloves, goggles, apron and mask
- Pour the required quantity of the ready-to-use solution into a tray so that the equipment to be disinfected is fully immersed

#### Use :

- 1- Immerse the medical device to be treated
- 2- Circulate the disinfectant inside the channels for endoscopes
- 3- Contact time: 10 minutes minimum
- 4- Rinse thoroughly with microbiologically controlled water

Renew the bath every 7 days maximum. Check the stability of the bath at least every morning before the first use with a PERATEST 500 strip (TDS 760\_P9H+).

### DISINFECTING PROPERTIES

Activity	Standards	Contact time
<b>Bactericidal</b>	EN 13727	1 minute
	EN 14561	5 minutes
<b>Fungicidal</b>	EN 13624 / EN 14562	5 minutes
<b>Virucidal</b>	EN 14476	10 minutes
	EN 17111	5 minutes
<b>Mycobactericidal</b>	EN 14348 / EN 14563	10 minutes
<b>Sporicidal</b>	EN 14347	10 minutes
	EN 17126	10 minutes
<b>Sporicidal against C. difficile</b>	EN 17126	10 minutes

### PRECAUTIONS OF USE

Refer to the Safety Data Sheet (SDS 010474)

### FEATURES

- Aspect : **Liquid**
- Colour : **Colourless**
- Pure pH : **2,5**
- Expiration date : **1 year**

### COMPOSITION

- Peracetic acid
- Hydrogen peroxide
- Acetic acid

### PACKAGING

- 1 Box of 4 x 5 L cans .....ref. **1047420**



Hazardous product - follow the precautions of use. Always read the label and product information before use. Class IIb medical device